



Anti-RPT2 (aa 1-100) polyclonal antibody (CABT-BL3195)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant (His-tagged) fusion protein: MGQGVSSGQD KKKKKGSNQG PKYEPPVQSK FGRKKRKGGP ATAELPNYIY PSTRCKLKLL RMERIKDHLL LEEEFVSNSE ILKPFEKKQE EEKKQLEEIR, corresponding to N terminal amino acids 1-100 of <i>S. cerevisiae</i> RPT2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	<i>S. cerevisiae</i>
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cytoplasmic and Nuclear
Format	Liquid
Buffer	Whole serum
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [851557](#)

Protein Refseq [NP_010277](#)

UniProt ID [P40327](#)
