



Anti-RPT2a (full length) polyclonal antibody (CABT-BL3196)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Full length protein produced in E.coli.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Arabidopsis thaliana
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cytoplasmic and Nuclear
Format	Liquid
Size	100 µl
Buffer	Whole serum
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	26S proteasome AAA-ATPase subunit RPT2a (RPT2a) mRNA
---------------------	--

GENE INFORMATION

Entrez Gene ID [829025](#)

Protein Refseq [NP_194633](#)

UniProt ID [Q9SZD4](#)
