



Anti-RPT3 (aa 1-100) polyclonal antibody (CABT-BL3197)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	His6 tagged recombinant fragment: MEELGIVTPVEKAVEEKPAVKSYASLLAQLNGTVNNSAL SNVNSDIYFKLKKLEKEYELLTLQEDYIKDEQRHLKRELK RAQEEVKRIQSVPLVIGQFL, corresponding to N terminal amino acids 1-100 of yeast, <i>S. cerevisiae</i> Rpt3 (YTA2) protein sequence.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	<i>S. cerevisiae</i>
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cytoplasmic and nuclear
Format	Liquid
Buffer	Whole serum.
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [852003](#)

Protein Refseq [NP_010682](#)

UniProt ID [P33298](#)
