



Anti-RRH (internal region) polyclonal antibody (CABT-BL3199)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide derived from an internal sequence of Human Peropsin
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Cellular Localization	Cell Membrane; Multi-pass membrane protein
Format	Liquid
Buffer	50% Glycerol, PBS, 150mM Sodium chloride, pH 7.4
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	Opsins are members of the guanine nucleotide-binding protein (G protein)-coupled receptor superfamily. This gene belongs to the seven-exon subfamily of mammalian opsin genes that includes opsin 5 and retinal G protein coupled receptor.
---------------------	---

GENE INFORMATION

Entrez Gene ID	10692
Protein Refseq	NP_006574
UniProt ID	Q14718