



Anti-RSPH10B (aa 288-337) polyclonal antibody (CABT-BL3201)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide, corresponding to a region within internal sequence amino acids 288-337 (FVNGYRHGRGKFYYASGAMYDGEWVSNKKHGMGRLTFKN GRVYEGAFSND) of Human RSPH10B
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	222967
Protein Refseq	NP_775836

UniProt ID

[B2RC85](#)
