



Anti-RSPH3 (aa 45-134) polyclonal antibody (CABT-BL3203)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to amino acids of human RSPH3.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	WB, IHC
Format	Liquid
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	RSPH3 (radial spoke 3 homolog (Chlamydomonas)) is a protein-coding gene. Diseases associated with RSPH3 include sleeping sickness, and cryptorchidism.
--------------	--

GENE INFORMATION

Entrez Gene ID [83861](#)

Protein Refseq [NP_114130.3](#)

UniProt ID [Q86UC2](#)
