



Anti-RTP2 (aa 130-160) polyclonal antibody (CABT-BL3210)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide conjugated to KLH, corresponding to a region within internal sequence amino acids 130-160 of Human RTP2
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cell membrane; Single-pass type III membrane protein. Note: Effective cell surface expression depends upon interaction with olfactory receptors
Format	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

GENE INFORMATION

Entrez Gene ID [344892](#)

Protein Refseq [NP_001004312](#)

UniProt ID [Q5QGT7](#)
