



Anti-RUSC1-AS1 (aa 143-192) polyclonal antibody (CABT-BL3219)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within internal sequence amino acids 143-192 (APLSCPAPRA QVHRSTPMGR ALLTRVLLEP LRPWACPRLP RSPPGGAQSG) of Human C1orf104.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	284618
Protein Refseq	NP_001034606

UniProt ID

[Q66K80](#)
