



Anti-SARA (aa 341-363) polyclonal antibody (CABT-BL3229)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide: DQQMPDDEKS CSASPVTEEEESVR conjugated to KLH, corresponding to Internal sequence amino acids 341-363 of Fruit fly (<i>Drosophila melanogaster</i>) SARA
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	<i>Drosophila</i>
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC, IHC-Fr, IP, WB
Cellular Localization	Cytoplasm. Early endosome membrane.
Format	Liquid
Buffer	Antibody stabilization buffer, with preservatives, consisting of Hepes, NaCl, KCl, glycerol, BSA and chaotropic agents. Preservative: 0.02% Sodium Azide
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [44263](#)

Protein Refseq [NP_524729](#)

UniProt ID [Q7K9H6](#)
