



# Anti-SAYSD1 polyclonal antibody (CABT-BL3237)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                    |                                                                                                     |
|--------------------|-----------------------------------------------------------------------------------------------------|
| Immunogen          | Recombinant protein corresponding to amino acids of human C6orf64.                                  |
| Isotype            | IgG                                                                                                 |
| Source/Host        | Rabbit                                                                                              |
| Species Reactivity | Human, Mouse, Rat                                                                                   |
| Purification       | Antigen affinity purification                                                                       |
| Conjugate          | Unconjugated                                                                                        |
| Applications       | WB, IHC                                                                                             |
| Format             | Liquid                                                                                              |
| Buffer             | In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)                                                   |
| Preservative       | 0.02% Sodium Azide                                                                                  |
| Storage            | Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing. |

## GENE INFORMATION

|                |                             |
|----------------|-----------------------------|
| Entrez Gene ID | <a href="#">55776</a>       |
| Protein Refseq | <a href="#">NP_060792.1</a> |
| UniProt ID     | <a href="#">Q9NPB0</a>      |