



Anti-SCAF1 (aa 1203-1216) polyclonal antibody (CABT-BL3239)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide conjugated to BSA, corresponding to amino acids 1203-1216 of Human SRA1
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Rat, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IP, ICC/IF
Cellular Localization	Nuclear
Format	Liquid
Buffer	50% Ammonium sulphate, PBS
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C (add glycerol to a final volume of 50% for extra stability). Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [58506](#)

Protein Refseq [NP_067051](#)

UniProt ID [Q9H7N4](#)
