



Anti-SCAMP5 (aa 217-233) polyclonal antibody (CABT-BL3243)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic Peptide: CNQPQTQYSATPNYTYSN. Immunizing peptide corresponds to amino acid residues 217-233 from rat SCAMP 5. This sequence is 94% conserved in mouse.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Rat, Cow
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ICC/IF
Cellular Localization	Cell membrane. Golgi apparatus membrane. Golgi apparatus; trans-Golgi network membrane. Recycling endosome membrane. Cytoplasmic vesicle; secretory vesicle; synaptic vesicle membrane. Mainly localizes in Golgi apparatus membrane. Upon calcium-triggered
Format	Liquid
Size	100 µg
Buffer	PBS, 1 mg/ml BSA and 0.05% sodium azide
Preservative	0.05% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

Secretory carrier-associated membrane protein 5 is a protein that in humans is encoded by the SCAMP5 gene.

GENE INFORMATION

Entrez Gene ID	56807
Protein Refseq	NP_064666
UniProt ID	Q9JKD3