



Anti-SCAPER (aa 1-50) polyclonal antibody (CABT-BL3248)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within N terminal amino acids 1-50 (MMASFQRSNS HDKVRRIVAE EGRTARNLIA WSVPLESKDD DGKPKCQTGG) of Human SCAPER.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Endoplasmic reticulum. Nucleus. Predominantly located in the endoplasmic reticulum, only a small portion is detected in the nucleus.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [49855](#)

Protein Refseq [NP_065894](#)

UniProt ID [Q9BY12](#)
