



Anti-SCD5 polyclonal antibody (CABT-BL3251)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to amino acids of human SCD5.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	IHC
Format	Liquid
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Fatty acid delta-9-desaturase that introduces a double bond in fatty acyl-coenzyme A at the delta-9 position.
---------------------	---

GENE INFORMATION

Entrez Gene ID [79966](#)

Protein Refseq [NP_079182.2](#)

UniProt ID [Q86SK9](#)
