



# Anti-SCD5 (middle region) polyclonal antibody (CABT-BL3252)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	This antibody is specific for the Middle Region of the target protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Buffer	20 mM Potassium Phosphate, 150 mM Sodium Chloride, pH 7.0
Preservative	None
Storage	Store at -20°C. Avoid freeze-thaw cycles.

## BACKGROUND

Introduction	Fatty acid delta-9-desaturase that introduces a double bond in fatty acyl-coenzyme A at the delta-9 position.
--------------	---

## GENE INFORMATION

Entrez Gene ID	<a href="#">79966</a>
----------------	-----------------------

Protein Refseq [NP\\_079182.2](#)

---

UniProt ID [Q86SK9](#)

---