



# Anti-Scn2a1 polyclonal antibody (CABT-BL3257)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Purified peptide of alpha subunit of rat Nav1.3 ( type III voltage-gated sodium channel VGSC)
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Rat
<b>Purification</b>	Affinity Purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	IHC, WB
<b>Preservative</b>	None
<b>Storage</b>	Maintain lyophilized material at -20°C for up to 23 months after date of receipt. After reconstitution maintain at -20°C in undiluted aliquots for up to 6 months. Avoid repeated freeze/thaw cycles.

## GENE INFORMATION

<b>Entrez Gene ID</b>	<a href="#">110876</a>
<b>Protein Refseq</b>	<a href="#">NP_001092768.1</a>
<b>UniProt ID</b>	<a href="#">B1AWN6</a>