



# Anti-SCRIN3 (internal region) polyclonal antibody (CABT-BL3260)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Synthetic peptide conjugated to KLH, corresponding to an internal region of Human SCRIN3
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Immunogen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, ELISA, IHC-P, ICC
<b>Format</b>	Liquid
<b>Buffer</b>	99% PBS
<b>Preservative</b>	0.02% Sodium Azide
<b>Storage</b>	Store at +4°C.

## BACKGROUND

<b>Introduction</b>	SCRIN3 (secernin 3) is a protein-coding gene. GO annotations related to this gene include dipeptidase activity. An important paralog of this gene is SCRIN2.
---------------------	--

## GENE INFORMATION

Entrez Gene ID [79634](#)

---

Protein Refseq [NP\\_078859](#)

---

UniProt ID [Q0VDG4](#)

---