



Anti-SCRIN3 polyclonal antibody (CABT-BL3261)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to amino acids of human SCRIN3.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	IHC
Format	Liquid
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	SCRIN3 (secernin 3) is a protein-coding gene. GO annotations related to this gene include dipeptidase activity. An important paralog of this gene is SCRIN2.
--------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------

GENE INFORMATION

Entrez Gene ID [79634](#)

Protein Refseq [NP_001180457](#)

UniProt ID [Q0VDG4](#)
