



Anti-SCRT2 (aa 258-307) polyclonal antibody (CABT-BL3262)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within internal sequence amino acids 258-307 (HMQTHSAFKH YRCRQCDKSF ALKSYLHKHC EAACAKAAEP PPPTPAGPAS) of Human Scratch 2
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Cellular Localization	Nuclear
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [85508](#)

Protein Refseq [NP_149120](#)

UniProt ID [Q9NQ03](#)
