



Anti-SEC14L3 (aa 271-320) polyclonal antibody (CABT-BL3270)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within C terminal amino acids 271-320 (RDQVKTQYEH SVQISRGS SH QVEYEILFPG CVLRWQFSSD GADIGFGVFL) of Rat Sec14l3 (NP_072130).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Rat
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	PBS, 2% Sucrose
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	64543
Protein Refseq	NP_072130

UniProt ID

[B1WC82](#)
