



Anti-SERHL2 (aa 237-267) polyclonal antibody (CABT-BL3275)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide conjugated to KLH, corresponding to a region within C terminal amino acids 237-267 of Human SERHL2
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cytoplasm, perinuclear region. Peroxisome. Note: Concentrated in perinuclear vesicles.
Format	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

GENE INFORMATION

Entrez Gene ID	253190
Protein Refseq	NP_055324

UniProt ID

[Q9H4I8](#)
