



Anti-SFPI1 (N-terminal) polyclonal antibody (CABT-BL3285)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to N-terminal residues of mouse PU.1/Spi1
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Cellular Localization	Nucleus.
Format	Liquid
Buffer	50% Glycerol, PBS
Preservative	0.01% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	20375
Protein Refseq	NP_035485

UniProt ID

[P17433](#)
