



Anti-SKD (aa 519-618) polyclonal antibody (CABT-BL3296)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Fusion protein: ADGDENEQNKPPQDSKQLVQQQIQQQQQRQKLWNFVDPMQ KAPCICKHLGNTPHGTTPHGGASTYSRNSLGGDSSMPVAS VESPATPAPSPHPNSAHSQP, corresponding to amino acids 519/618 of Drosophila melanogaster Skuld.
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Fruit fly
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nuclear
Format	Liquid
Buffer	50% Glycerol
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [43906](#)

Protein Refseq [NP_730590](#)

UniProt ID [E6PBY1](#)
