



Anti-SKIV2L2 polyclonal antibody (CABT-BL3297)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide from a region within amino acids 100-150 of Human SKIV2L2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Cellular Localization	Nucleus, nucleolus.
Format	Liquid
Buffer	0.1% BSA, Tris buffered saline
Preservative	0.09% Sodium Azide
Storage	Store at +4°C.

BACKGROUND

Introduction	Superkiller viralicidic activity 2-like 2 is a protein that in humans is encoded by the SKIV2L2 gene.
---------------------	---

GENE INFORMATION

Entrez Gene ID [23517](#)

Protein Refseq [NP_056175](#)

UniProt ID [P42285](#)
