



Anti-SLC10A3 (aa 287-336) polyclonal antibody (CABT-BL3299)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within internal sequence amino acids 287-336 (MTFLSTVAAT GFLPLSSAIY SRLLSIHETL HVPISKILGT LLFIAPIAV) of Human SLC10A3.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Membrane.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	This gene maps to a GC-rich region of the X chromosome and was identified by its proximity to a CpG island. It is thought to be a housekeeping gene. Alternative splicing results in multiple transcript variants encoding distinct isoforms.
---------------------	---

GENE INFORMATION

Entrez Gene ID [8273](#)

Protein Refseq [NP_001135864](#)

UniProt ID [P09131](#)
