



Anti-SLC10A5 (aa 361-410) polyclonal antibody (CABT-BL3300)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A region within synthetic peptide: PLPVCKTVAI ESGMLNSFLA LAVIQLSFPQ SKANLASVAP FTVAMCSGCE, corresponding to C terminal amino acids 361-410 of Human SLC10A5
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA, WB, IHC-P
Cellular Localization	Membrane; Multi-pass membrane protein
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	347051
Protein Refseq	NP_001010893

UniProt ID

[Q5PT55](#)
