



Anti-SLC22A10 (aa 431-480) polyclonal antibody (CABT-BL3307)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within internal sequence amino acids 431-480 (QTLRVALACL GIGCSAATFS SVAVHFIELI PTVLRARASG IDLTASRIGA) of Human SLC22A10.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Membrane.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [387775](#)

Protein Refseq [NP_001034841](#)

UniProt ID [Q63ZE4](#)
