



Anti-SLC25A44 (aa 264-313) polyclonal antibody (CABT-BL3317)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within C terminal amino acids 264-313 (QLMAEEGPWG LMKGLSARII SATPSTIVIV VGYESLKKLS LRPELVDSRH) of Human SLC25A44.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Mitochondrion inner membrane.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	SLC25A44 belongs to the SLC25 family of mitochondrial carrier proteins (Haitina et al., 2006 (PubMed;16949250)).
---------------------	--

GENE INFORMATION

Entrez Gene ID [9673](#)

Protein Refseq [NP_055470](#)

UniProt ID [E9PGQ0](#)
