



Anti-SLC26A11 polyclonal antibody (CABT-BL3321)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to amino acids of human SLC26A11.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	WB, IHC
Format	Liquid
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	The members of solute carrier family 26 (SLC26) are multifunctional anion exchangers with diverse substrate specificity. They are mainly expressed on epithelial cells and play a central role in anion secretion and absorption. SLC26A11 is widely distribute
---------------------	---

GENE INFORMATION

Entrez Gene ID [284129](#)

Protein Refseq [NP_001159819](#)

UniProt ID [Q86WA9](#)
