



Anti-SLC30A2 (internal region) polyclonal antibody (CABT-BL3325)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide derived from an internal sequence of Human SLC30A2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Cellular Localization	Lysosome membrane and Vacuole membrane.
Format	Liquid
Buffer	50% Glycerol, 0.88% Sodium chloride, 49.1% PBS pH: 7.40 Note: PBS is without Mg2+ and Ca2+.
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Stable for 12 months at -20°C

GENE INFORMATION

Entrez Gene ID	7780
Protein Refseq	NP_001004434

UniProt ID

[B3KSN7](#)
