



# Anti-SLC36A1 polyclonal antibody (CABT-BL3341)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to amino acids of human SLC36A1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	IHC
Format	Liquid
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

Introduction	This gene encodes a member of the eukaryote-specific amino acid/auxin permease (AAP) 1 transporter family. The encoded protein functions as a proton-dependent, small amino acid transporter. This gene is clustered with related family members on chromosome 5q33.1.
--------------	--

# GENE INFORMATION

Entrez Gene ID	<a href="#">206358</a>
Protein Refseq	<a href="#">NP_510968</a>
UniProt ID	<a href="#">Q7Z2H8</a>