



Anti-SLC41A1 (N-terminal) polyclonal antibody (CABT-BL3346)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to a region within the N-terminal sequence 35-84 (TSEFLGPDGAGVEVVIESRANAKGVREEDALLENGSQSN ESDDVSTDRGP) of Human SLC41A1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Cellular Localization	Cell membrane; Multi-pass membrane protein.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [419837](#)

UniProt ID

[E1C334](#)
