



Anti-SLC44A3 (aa 287-336) polyclonal antibody (CABT-BL3349)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to internal sequence amino acids 287-336 (TFAILIFFWV LWVAVLLSLG TAGAAQVM EG GQVEYKPLSG IRYMWSYHLI) of Human SLC44A3
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Reconstitution	Add 50ul distilled water for a final concentration of 1 mg/ml.
Cellular Localization	Membrane
Format	Lyophilised
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [126969](#)

Protein Refseq [NP_001107578](#)

UniProt ID [Q8N4M1](#)
