



# Anti-SLC46A3 polyclonal antibody (CABT-BL3353)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to amino acids of human SLC46A3.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	WB, IHC
Format	Liquid
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

Introduction	Folate, a member of the vitamin B complex, is obtained solely from dietary sources. It is taken up from epithelial cells into tissues such as the kidney, liver, and intestine via folate transport systems. Members of solute carrier family 46 (SLC46) are fo
--------------	---

# GENE INFORMATION

Entrez Gene ID [283537](#)

---

Protein Refseq [NP\\_001129391](#)

---

UniProt ID [Q7Z3Q1](#)

---