



Anti-SLC7A6 (C-terminal) polyclonal antibody (CABT-BL3360)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide selected from the C-terminal region of Human SLC7A6 conjugated to KLH.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Basolateral cell membrane; Multi-pass membrane protein.
Format	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

GENE INFORMATION

Entrez Gene ID	9057
Protein Refseq	NP_001070253

UniProt ID

[A0A024R719](#)
