



Anti-SLC9A2 polyclonal antibody (CABT-BL3363)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to amino acids of human SLC9A2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	IHC
Format	Liquid
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction Sodium–hydrogen exchanger 2 is a protein that in humans is encoded by the SLC9A2 gene.

GENE INFORMATION

Entrez Gene ID [6549](#)

Protein Refseq [NP_003039](#)

UniProt ID [Q9UBY0](#)
