



Anti-SLC9B1 (aa 1-50) polyclonal antibody (CABT-BL3366)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within N terminal amino acids 1-50 (MHTTESKNEH LEDENFQTST TPQSLIDPNN TAHEETKTVL SDTEEIKPQT) of Human NHEDC1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Membrane.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	The protein encoded by this gene is a sodium/hydrogen exchanger and transmembrane protein. Highly conserved orthologs of this gene have been found in other mammalian species.
---------------------	--

The expression of this gene may be limited to testis. Multiple transcript variants encoding different isoforms have been found for this gene.

GENE INFORMATION

Entrez Gene ID	150159
----------------	------------------------

Protein Refseq	NP_631912
----------------	---------------------------

UniProt ID	Q4ZJI4
------------	------------------------
