



Anti-SLC01A2 (aa 451-500) polyclonal antibody (CABT-BL3368)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within internal sequence amino acids 451-500 (VCGNNGLSYL SACLAGECTS IGTGINMVFQ NCSCIQTSGN SSAVLGLCDK) of Human SLC01A2 (AAD52694).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cell membrane.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	This gene encodes a sodium-independent transporter which mediates cellular uptake of organic ions in the liver. Its substrates include bile acids, bromosulphophthalein, and some steroidal
---------------------	---

compounds. The protein is a member of the SLC21A family of solute carriers. Alternative splicing of this gene results in multiple transcript variants.

GENE INFORMATION

Entrez Gene ID	6579
----------------	----------------------

Protein Refseq	NP_066580
----------------	---------------------------

UniProt ID	A0A024RAT5
------------	----------------------------
