



Anti-SLCO1B7 (aa 576-625) polyclonal antibody (CABT-BL3369)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide within residues: RAFFGLKVAL IFPVLVLLTV FIFVVRKKSH GKDTKVLENE RQVMDEANLE, corresponding to internal sequence amino acids 576-625 of Human LST-3TM12.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	338821
Protein Refseq	NP_001009562

UniProt ID

[G3V0H7](#)
