



# Anti-SLCO2B1 polyclonal antibody (CABT-BL3370)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Recombinant protein corresponding to amino acids of human SLCO2B1.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Antigen affinity purification
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, IHC
<b>Format</b>	Liquid
<b>Buffer</b>	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
<b>Preservative</b>	0.02% Sodium Azide
<b>Storage</b>	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

<b>Introduction</b>	This locus encodes a member of the organic anion-transporting polypeptide family of membrane proteins. The protein encoded by this locus may function in regulation of placental uptake of sulfated steroids. Alternatively spliced transcript variants have been described.
---------------------	--

# GENE INFORMATION

Entrez Gene ID [11309](#)

---

Protein Refseq [NP\\_001138683](#)

---

UniProt ID [O94956](#)

---