



Anti-SLFN2 (C-terminal) polyclonal antibody (CABT-BL3373)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to C-terminal residues of mouse SLFN2 or to N-terminal residues of mouse SLFN3 or to N-terminal residues of mouse SLFN4
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Format	Liquid
Buffer	50% Glycerol, PBS
Preservative	0.01% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	20556
Protein Refseq	NP_035538
UniProt ID	Q9Z0I6