



Anti-SMG9 (aa 75-125) polyclonal antibody (CABT-BL3385)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide designed within residues corresponding to amino acids 75 - 125 of Human C19orf61.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IP, WB
Format	Liquid
Size	100 µl
Buffer	Tris citrate/phosphate buffer, pH 7-8
Preservative	0.09% Sodium Azide
Storage	Store at +4°C.

GENE INFORMATION

Entrez Gene ID	56006
Protein Refseq	NP_061981

UniProt ID

[A0A024R0Q0](#)
