



Anti-SNAPC5 (aa 10-59) polyclonal antibody (CABT-BL3389)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within N terminal amino acids 10-59 (KEEETLLRLK AALHDQLNRL KVEELALQSM ISSRRGDEML SSHTVPEQSH) of Human SNAPC5
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nuclear
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	snRNA-activating protein complex subunit 5 is a protein that in humans is encoded by the SNAPC5 gene.
---------------------	---

GENE INFORMATION

Entrez Gene ID [10302](#)

Protein Refseq [NP_006040](#)

UniProt ID [A0A024R5X1](#)
