



Anti-SNF2 (aa 2-32) polyclonal antibody (CABT-BL3393)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Fusion protein: NIPQRQFSNEEVNRCYLRLWQHRLNEHGMNAP, corresponding to N terminal amino acids 2/32 of <i>Saccharomyces cerevisiae</i> SWI2 / SNF2
Isotype	IgG
Source/Host	Mouse
Species Reactivity	<i>S. cerevisiae</i>
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nuclear
Format	Liquid
Buffer	50% Glycerol
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [854465](#)

Protein Refseq [NP_014933](#)

UniProt ID [P22082](#)
