



Anti-SOBP polyclonal antibody (CABT-BL3403)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to amino acids of human SOBP.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	WB, IHC
Format	Liquid
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	The protein encoded by this gene is a nuclear zinc finger protein that is involved in development of the cochlea. Defects in this gene have also been linked to intellectual disability.
---------------------	--

GENE INFORMATION

Entrez Gene ID [55084](#)

Protein Refseq [NP_060483](#)

UniProt ID [A7XYQ1](#)
