



Anti-SOGA2 polyclonal antibody (CABT-BL3404)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, WB, IHC
Format	Lyophilized Powder
Buffer	PBS with 0.02% Sodium Azide, 1%BSA,PH7.2
Preservative	0.02% Sodium Azide
Storage	-20°C, Avoid freeze / thaw cycles

GENE INFORMATION

Entrez Gene ID	23255
Protein Refseq	NP_056025.2
UniProt ID	Q9Y4B5