



Anti-SOHLH1 (aa 1-100) polyclonal antibody (CABT-BL3405)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide conjugated to KLH derived from within residues 1 - 100 of Mouse Sohlh1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, WB
Cellular Localization	Cytoplasm. Nucleus.
Format	Liquid
Buffer	1% BSA, PBS, pH 7.4
Preservative	0.02% Sodium Azide
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Probable transcription factor required during spermatogenesis and oogenesis.
---------------------	--

GENE INFORMATION

Entrez Gene ID [227631](#)

Protein Refseq [NP_001001714](#)

UniProt ID [Q6IUP1](#)
