



Anti-SPATA24 (aa 151-200) polyclonal antibody (CABT-BL3421)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within C terminal amino acids 151-200 (LQQVISQQKQ IFRNHMSDFR IQKQQESYMA QVLDQKHKKA SGTRQARSHQ) of Human SPATA24.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cytoplasm. Nucleus; nucleolus. Nucleus; nucleoplasm. Associated with chromatin.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [202051](#)

Protein Refseq [NP_919272](#)

UniProt ID [Q86W54](#)
