



Anti-SPCC74.02c (aa 454-553) polyclonal antibody (CABT-BL3422)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Fusion protein: VSNNREELPSFRKRSSSRQSEEIASLQAENAVFSSDPAS NDEKTGNKKRKKKSVSWKPDNDLVQVKFIESLNEEGAASV KTPHIYGNARDMDRQEARVA, corresponding to amino acids 454/553 of <i>S. pombe</i> PAXH.
Isotype	IgG
Source/Host	Mouse
Species Reactivity	<i>S. pombe</i>
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	50% Glycerol
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	2539414
Protein Refseq	NP_588375

UniProt ID

[O74535](#)
