



Anti-SPICE1 (aa 36-85) polyclonal antibody (CABT-BL3428)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within N terminal amino acids 36-85 (VTDLTVHRAT PEDLVRRHEI HKSKNRALVH WELQEKALKR KWRKQKPETL) of human CCDC52.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Reconstitution	Add 50 µl of distilled water for a final concentration of 1 mg/ml.
Cellular Localization	Cytoplasm; cytoskeleton; centrosome; centriole. Cytoplasm; cytoskeleton; spindle.
Format	Lyophilised
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [152185](#)

Protein Refseq [NP_653319](#)

UniProt ID [Q8N0Z3](#)
